



PHILIPS

IVUS



Coronary IVUS

Image Guided Therapy – Devices
2021

innovation ✨ you



PHILIPS

IVUS

Contents

Philips IGT-D	03
What is CAD	09
What is IVUS	13
Why perform IVUS in CAD	20
Our products	44
Clinical data	54



Philips IGT-D

Image Guided Therapy - Devices



Ready to take on the healthcare challenge?



At Philips, we take a **holistic view of people's health journeys**, starting with healthy living and prevention, precision diagnosis and personalized treatment, through to care in the home – where the cycle to healthy living begins again.



Seamless care – designed to meet our customers' pressing needs



Simplifying
data and insight gathering



Removing
excess costs



Driving
improved treatment
and outcomes



Giving
patients and staff a
better experience



Our vision for Image Guided Therapy

We're leading in thought and action so that our patients can get home faster and healthier, and stay there, **living the lives they love.**



IGT-devices

An extensive and differentiated portfolio for a new value proposition

Lead management & EP

- Lead locking devices
- Mechanical and laser extraction sheaths
- Rescue occlusion balloon

Superficial femoral & popliteal arteries

- PVIVUS
- Laser atherectomy Atherectomy
- Support catheters
- Drug coated balloons
- Scoring balloons

Coronary

- Coronary IVUS
- Pressure & flowwires
- Scoring balloons
- Laser atherectomy

Dialysis access


- Peripheral vascular IVUS
- Scoring balloon

Venous

- Peripheral vascular IVUS

Below the knee

- Laser atherectomy, Atherectomy
- Peripheral vascular IVUS
- Support catheters



— Coronary — Peripheral — Lead management

Image Guided Therapy

Coronary devices

IVUS imaging	Physiology	Therapies	Advanced imaging solutions
<p>Eagle Eye Platinum, Refinity & Revolution</p>	<p>Verrata Plus & Verrata</p>	<p>ELCA, AngioSculpt & AngioSculptX</p>	<p>SyncVision</p>

Your partner in complex PCI



What is CAD



Coronary arteries


There are 2 main coronary arteries of the heart

The left main coronary artery (LMCA) supplies blood to the cardiac muscle on the left side of the heart (the left ventricle and left atrium). The LMCA divides into branches:

- The left anterior descending artery branches off the left main coronary artery and supplies blood to the front of the left side of the heart.
- The circumflex artery also branches off the left main coronary artery and encircles the heart muscle. This artery supplies blood to the outer side and back of the heart.

The right coronary artery (RCA) supplies blood to the cardiac muscle on the right side of the heart (right ventricle and the right atrium) and the SA (sinoatrial) and AV (atrioventricular) nodes, which regulate the heart rhythm. The right coronary artery divides into smaller branches, including the right posterior descending artery and the acute marginal artery. Together with the left anterior descending artery, the right coronary artery helps supply blood to the middle or septum of the heart.

- Smaller branches of the coronary arteries include: obtuse marginal (OM), septal perforator (SP), and diagonals.






Structure of an artery

The wall of an artery consists of three layers

The innermost layer (tunica intima), is simple endothelium surrounded by a connective tissue basement membrane with elastic fibers.

The middle layer (tunica media), is primarily smooth muscle and is quite often the thickest layer, provides support for the vessel but also changes vessel diameter to regulate blood flow and blood pressure.

The outermost layer (tunica externa or tunica adventitia), attaches the vessel to the surrounding tissue. This layer is connective tissue with varying amounts of elastic and collagenous fibers. The connective tissue in this layer is quite dense where it is adjacent to the Tunica media, but it changes to loose connective tissue near the periphery of the vessel.





Coronary artery disease (CAD)


What is CAD?

Coronary artery disease develops when the major blood vessels that supply your heart with blood, oxygen and nutrients (coronary arteries) become damaged or diseased. Cholesterol-containing deposits (plaque) in your arteries and inflammation are usually to blame for coronary artery disease.

When plaque builds up, it narrows your coronary arteries, decreasing blood flow to your heart. Eventually, the decreased blood flow may cause chest pain (angina), shortness of breath, or other coronary artery disease signs and symptoms. A complete blockage can cause a heart attack.

Coronary artery disease is thought to begin with damage or injury to the inner layer of a coronary artery, sometimes as early as childhood. The damage may be caused by various factors, including:

- Smoking
- High blood pressure
- High cholesterol
- Diabetes or insulin resistance
- Sedentary lifestyle







What is IVUS

What is IVUS?

Intravascular ultrasound

- When electricity stimulates a crystal, sound waves are released
- The sound hits tissue and is reflected to a varying degree
- Crystals then receive those reflected signals and code them back into electrical impulses
- Those impulses can construct a sophisticated image using computing power



Adventitia
Whitering

Intima

Media
Blackring

Dead Zone
Catheter

Two ways to acquire IVUS signals


Digital IVUS - Eagle Eye Platinum

A series of fixed 64 ultrasound devices that scan sections of the vessel in turn, whilst computers to build the image, 20MHz greyscale.




Rotational IVUS - Refinity & Revolution

A single IVUS catheter that spins within a sheath and uses computers to build the image in real-time from the sections that come in bit by bit, 45MHz greyscale.





MHz – megahertz

MHz refers to the frequency of the soundwave. A higher frequency will increase the resolution but reduce the depth penetration.



Plaque geometry

Atherosclerotic plaques typically have different geometries. The two types of geometry you will see are:



Eccentric plaques are distributed non-circumferentially in the vessel; this makes the assessment of disease by angiography especially prone to underestimation or overestimation depending on the angle of view.


Concentric plaques are distributed circumferentially in the vessel and tend to occur in areas of negative remodeling; use of angiography alone could result in too large a stent diameter.



Four different types of plaque


1. Soft (fatty)

Echolucent light gray flecks




2. Fibrous

Echogenic, light gray with white surfaces




3. Calcified

Highly echogenic, white areas with shade





4. Mixed plaque

Mixed plaque is a combination of tissues of varying echogenicity



Two types of IVUS image

Rotational & digital





Grayscale levels

Plaque composition

Calcified

Fibrotic

Fibrous

Soft



Calcium is indicated by very bright areas with acoustic shadowing that blocks out the image behind. This shadow occurs because the high density of calcium dampens the ultrasound echo.

Fibrous plaques have an intermediate echogenicity between soft (echolucent) atheroma and highly echogenic calcified plaques. Fibrous plaques exhibit little or no acoustic shadowing.

Plaque definition should be determined by reviewing the entire segment, not just the image slice with the smallest lumen CSA.



Why perform IVUS in CAD



What beneficial information does intravascular ultrasound provide?

- Vessel sizing
- Length of lesions
- Type of plaque
- Measurement of plaque/stenosis
- Optimized treatment strategies based on morphology and geometry
- Detection of dissections & complications with PCI





IVUS adds a third dimension

An angiogram only provides a 2D image of a vessel lumen, referred to as **lumen imaging**¹


- The clinical problem involves a 3D vessel and lesion

With 2D coronary angiography:

- Depending on the angle of the x-ray arm producing the angiogram, a value for % stenosis may be established that does not reflect the true narrowing of the vessel²
- An under- or over-estimation of a lesion's severity may result³

As a 3D modality, IVUS can accurately identify disease, assess its significance and provide valuable additional information:^{4,5}

- Assess plaque burden and lesion morphology^{4,6}
- Detect acute complications^{6,7}
- Provide accurate measurements for device selection⁴
- Assess adequacy of stent deployment⁵




1. McDaniel MC, et al. JACC Cardiovasc Interv. 2011;4(11):1155-67.
 2. Topol EJ, Nissen SE. Circulation. 1995;92(8):2333-42.
 3. Pijls NH, et al. J Am Coll Cardiol. 2010;56(3):177-84.
 4. McDaniel MC, et al. JACC Cardiovasc Interv. 2011;4(11):1155-67.

5. Bonello L, et al. Arch Cardiovasc Dis. 2009;102(2):143-51.
 6. Porto I, et al. Cardiovasc Ultrasound. 2004;2:18.
 7. Yip A, S J. Cardiovasc Diagn Ther. 2015;5(1):37-48.




Angiography vs IVUS






Angiography alone is not enough:

Ambiguous angiogram




- Dissection?
- Calcium?
- Thrombus?
- Stenosis?


Dissection



Calcium



Thrombus



Stenosis




Without IVUS, how would you know which condition to treat from this angiogram?



IVUS provides information to guide treatment decisions throughout the procedure


Pre-procedure

Identify disease location, assess significance and detect acute complications




Peri-procedure

Select appropriate treatment strategy dependent on lesion morphology and select appropriate devices



Post-procedure

Post-treatment assessment to optimize procedural outcomes



Pre-procedure IVUS

Plan the intervention

1. Identify disease & assess significance^{1,2}

Lesion severity

- Cross sectional lumen area
- Vessel diameter
- Plaque burden
- Calcification
- Sub-intimal thickening
- Thrombus

In case of ambiguous angiogram



- Diffuse reference vessel disease
- Lesion foreshortening
- Angulations
- Eccentricity
- Hazy lesions

Side-branches

- Presence of vascular remodeling

2. Detect acute complications^{2,3}

- Plaque rupture
- Spontaneous dissection
- Thrombus



1. McDaniel MC, et al. JACC Cardiovasc Interv. 2011;4(11):1155-67.
2. Porto I, et al. Cardiovasc Ultrasound. 2004;2:18.
3. Yip A, S J. Cardiovasc Diagn Ther. 2015;5(1):37-48.




Thrombus

Blood clot within the vessel

One of the hardest things to determine is thrombus. Many times it is missed altogether on angiography and can only be seen with IVUS. Thrombus will typically appear as an irregular shape within the vessel and may be moving with flow through the vessel.

It can be much easier to see in a moving IVUS video, rapid review, or by using ChromaFlo imaging than in a single still frame.

Stent thrombosis (ST) is a thrombotic occlusion of a coronary stent. This is usually an acute process in contrast to restenosis, which is a gradual narrowing of the stent lumen due to neointimal proliferation. Stent thrombosis often results in an acute coronary syndrome, while restenosis often results in symptoms of angina.





Thrombus

In-stent restenosis (ISR)

Plaque within stent

ISR can be classified as >50% diameter stenosis at stent's inside or edges, with different ISR morphologies predicting different TLR rates at 1-year^{1,2}. Several factors with both biological and/ or procedural causes contribute to the nature of ISR³



1st stent

2nd stent

Biologic causes

- Reaction to metal or polymer
- Drug resistance
- Thrombosis

Procedural causes

- Stent under-expansion/mal-apposition
- Stent fracture
- Edge trauma
- Geographical miss

1. Mehran R, Dangas G, Abizaid AS, et al. Angiographic patterns of in-stent restenosis: classification and implications for long-term outcome. *Circulation* 1999;100:1872–8.
 2. Solinas E, Dangas G, Kirtane AJ, et al. Angiographic patterns of drug eluting stent restenosis and one-year outcomes after treatment with repeated percutaneous coronary intervention. *Am J Cardiol* 2008;102:311–5.
 3. Mintz GS. "Clinical Utility of Intravascular Imaging and Physiology in Coronary Artery Disease". *J Am Coll Cardiol* 2014;64:207–22.


Dissection

Intimal tear resulting in the formation of a 'false' and 'true' lumen

Tears in the plaque that are parallel to the vessel wall with visualization of blood flow in the false lumen.



Can be difficult to detect live. Most often occurs at a juncture between varying densities:

- Stent/soft plaque
- Calcium to soft plaque
- Soft plaque to wall



Example of ChromaFlo imaging highlighting a dissection

Deep wall dissection



- Concentric lesion
- Percutaneous transluminal coronary angioplasty
- Shearing force opposite the calcified segment

Peri-procedure

Select the appropriate treatment strategy and use of devices



Assess lesion morphology and identify the need for pre-treatment¹

- Superficial or deep calcification
- Lesion tissue types to identify thin-cap fibroatheroma


Perform measurements to help select appropriate device(s)¹

- Vessel & lumen diameter
- Area & lesion length
- Lesion preparation
- Define appropriate landing zones
- Stent selection and sizing


IVUS allows you to select appropriate proximal and distal reference segments. Having a **defined, optimal PCI criteria** from these references, **significantly reduces MACE²** compared to angiography alone.




Angiography does not reveal vessel remodeling



Proximal reference



Lesion



Distal reference

1. McDaniel MC, et al. JACC Cardiovasc Interv. 2011;4(11):1155-67.
 2. Zhang J et al. The ULTIMATE trial. Intravascular Ultrasound-Guided Versus Angiography-Guided Implantation of Drug-Eluting Stent in All-Comers. Journal of the American College of Cardiology (2018), doi.org/10.1016/j.jacc.2018.09.013


Peri-procedure

Vessel measurements

Vessel diameter


Measure adventitia to adventitia (EEM to EEM)

- Maximum diameter
- Minimum diameter



Vessel/external elastic membrane cross-sectional area (EEM CSA)

- Vessel area measured at the site of external elastic membrane (outer boundary of media)
- Cross-sectional area inside of adventitia



Minimum diameter



Maximum diameter

EEM CSA



Peri-procedure

Lumen measurements



Lumendiameter

- Measured intima to intima
- Maximum diameter
- Minimum diameter

Minimal lumen diameter (MLD)

- Smallest lumen diameter within lesion segment

Lumen area

- Cross-sectional area inside of lumen

Minimal lumen area (MLA)


- Smallest Lumen CSA within lesion segment

Peri-procedure

Calculation of reference lumen area

Reference lumen area (RLA)

RLA is the average of proximal lumen CSA & distal lumen CSA





Peri-procedure

Calculation of % stenosis


Reference lumen area (RLA)

Average of proximal lumen CSA & distal lumen CSA



Minimal lumen area (MLA)

Smallest lumen CSA within lesion segment



$$\% \text{ stenosis} = \left[\frac{\text{RLA} - \text{MLA}}{\text{RLA}} \right] \times 100$$

Peri-procedure

Calculation of plaque burden

Remember


- Atheroma area (or plaque + media) H (EEM CSA – MLA)
- Plaque burden H $[(\text{Atheroma CSA}) / \text{EEM CSA}] \times 100$

Example*:

$$[(10.17 - 2.92) / 10.17] \times 100 =$$

71.29% plaque burden

% of plaque burden can be measured and calculated from a single IVUS frame.



EEM
Lumen

*Numbers have been rounded to one decimal place on IVUS image




Peri-procedure

PCI stent sizing

Several methods and techniques are used including:

- Stent sizing for remodelled vessels: "Mid-wall calculation" Using distal reference (vessel diameter + lumen diameter)/2 = approximate stent size.
- 80% of average (distal & proximal) reference diameters

Positively remodelled vessel





Post-procedure IVUS


Optimize procedural outcomes

Assess adequacy of stent expansion, apposition, mechanical complications and mechanisms:


- Geographical miss
- Major edge dissections
- Apposition
- Plaque protrusion




Undersized stent




Still undersized post high-pressure inflations



Minimum stent area above a certain threshold



Avoiding major malapposition



Absence of significant edge dissection




Post-procedure IVUS

Geographical miss (GM)

Inadequate stent deployment or suboptimal placement

Axial GM: Under or over inflation of a stent or balloon in which the ratio between the size of the largest balloon (at its maximum inflation pressure) and the reference vessel diameter was ≤ 0.9 or > 1.3 .

Longitudinal GM: A deployed stent which does not cover an injured or significantly diseased segment at one or both of the edges.



STLLR study reported: 66.5% of stents had “geographic miss,” associated with 3 times the MI and 2 times the TVR.¹

Stent optimization (STOP) study reported: IVUS-guidance increased the proportion of cases where optimal stent CSA (cross-sectional area) was achieved (21% & 54% vs 81%) following high pressure deployment and post dilation²

1. Costa et al. The Impact of Stent Deployment Techniques on Clinical Outcomes of Patients Treated With the CYPHER® Stent (S.T.L.L.R.). Am J Cardiol 2008 Jun 15;101(12):1704-11.

2. Rana O et al. The impact of routine intravascular ultrasound guided high-pressure post-dilation after drug-eluting stent deployment. J Invasive Cardiol 2014;26(12):640-646.





Post-procedure IVUS

Malapposed stents

Stent malapposition (SM): also referred to as incomplete stent apposition, is defined by the separation of at least one stent strut from the intimal surface of the arterial wall with evidence of blood behind the strut, without involvement of side branches.

This phenomenon, commonly identified by intravascular ultrasound (IVUS) imaging studies, may be detected early, at the time of stent implantation (and classified as acute) or later, at follow-up (therefore classified as late).¹




1. Karalis I, Ahmed TAHN, Jukema JW. Late acquired stent malapposition: why, when and how to handle? Heart 2012;98:1529-1536.

Post-procedure IVUS

Stent malapposition

It is possible to see malapposition with grayscale IVUS but ChromaFlo imaging enhances the differential between vessel and stent by highlighting the blood flow red.




Blood flow behind stent


Stent in contact with vessel wall

Mal-apposed stent


Use ChromaFlo for easy image interpretation




Mal-apposed stent



Well-apposed stent



Bifurcation
– absence of plaque in carina can influence provisional stenting strategy





Thrombus



Post-procedure IVUS

Well apposed stents

Both examples below show the stent well apposed against the vessel wall. No blood flow is observed between the stent struts and the intima.






Post-procedure IVUS

Opposed but undersized/under deployed stent

An example of a placed stent which appears to be well apposed, no space or blood flow is visible between the stent and the vessel wall but in this image the stent may be slightly undersized or under deployed.

It is only possible to determine if the stent was undersized or under deployed if you have all of the information pertaining to the placed stent.



Undersized or under deployed stent

Apposed stent

Post-procedure IVUS

Under-expanded stent



Looking at these two angiograms, it is hard to determine how well the stent inside the vessel has been deployed.

By using IVUS, displayed in the corresponding images, it is now possible to see that the angio on the top has a malapposed, under-expanded stent.

Provided with this information, the physician case went back in, post dilated the stent and ended up with the result you see in the bottom image.



It is extremely likely that if IVUS had not been performed, there is a good chance that the under-expansion of the stent may have been missed.

Lumen at proximal stent site with normal appearance



1: 2.9mm
2: 3.0mm

Final angio with optimal expansion



1: 3.5mm
2: 3.5mm



Our products

CORE precision guided therapy system

The choice of imaging and physiology on a single platform

CORE integrated & CORE mobile system

- Choice of multiple imaging and physiology tools on a single platform, only Philips offers digital IVUS and iFR
- Convenience of an integrated or mobile solution
- Intuitive interface for the ultimate in ease of use





IntraSight

Interventional applications platform

Smart. Simple. Seamless.

As the number of cath lab patients grows, so does the need to work smarter and faster.

Philips IntraSight offers you a comprehensive suite of clinically proven imaging, physiology and co-registration* tools on a modern, secure platform that will help you simplify complex interventions, speed routine procedures and improve lab efficiencies.

These best-in-class interventional tools ultimately allow you to go beyond the angiogram and complete your view of the target vessel so you can make fast, informed clinical decisions.

IntraSight is built on a modern, scalable platform that will be ready to provide you with new innovations and tools as they become available in the future.

For more information, go to www.philips.com/IntraSight

*Co-registration tools available within IntraSight 7 configuration via SyncVision




Coronary IVUS

Intravascular ultrasound

Rotational IVUS:


Refinity, Revolution & SpinVision



Quick and easy to use measurement, annotation and labelling tools


Digital IVUS:

Eagle Eye Platinum & Eagle Eye Platinum ST



ChromaFlo imaging:


Provides clear confirmation of stent apposition and lumensize






Refinity

Our next generation rotational IVUS catheter




- Low profile, 5Fr guide compatible
- Suitable for radial access
- GlyDx hydrophilic coating
- Soft and flexible tip with short distance transducer
- Smooth proximal to mid-shaft transition

Exceptional deliverability



Refinity ST insertion and recovery forces are significantly lower ($p < 0.001$) than OptiCross*

45 MHz, high resolution imaging




* Data on file at Philips. Bench testing was conducted with three rotational IVUS catheters: Refinity ST, Revolution, OptiCross through a tortuous model




Refinity & Revolution

Rotational IVUS catheters



Refinity



Revolution



Eagle Eye Platinum


Phased-array, digital IVUS catheter

- Fast, plug-and-play simplicity
- Quick, convenient length estimation
- ChromaFlo imaging
- Unique SyncVision IVUS & iFR co-registration compatibility
- Standard & short-tip variations




Eagle Eye Platinum

Phased-array, digital IVUS catheter



Eagle Eye Platinum




Eagle Eye Platinum ST

SyncVision advanced imaging solution

iFR & IVUS co-registration with angiographic enhancement

Advanced angiography: Angio+



- Vessel enhancement for clarity visualizing vessels
- QCA in real-time
- Device detection (stabilization & enhancement) to enable precise delivery

Advanced physiologic Imaging: iFR co-registration




- iFR drop is displayed on angiogram to identify regions of ischemia
- Length measurement without pullback device to size stent
- Estimate post-procedure iFR

Advanced IVUS imaging: IVUS co-registration



- Localization of IVUS with angiography to help avoid geographic miss
- Easy length/area/diameter measurements with manual pullback to size stent



“IVUS & iFR co-Registration help, with great ease of use, to address the PCI stent challenges of geographical mismatch, lesion location and PCI decision making.”

Prof. Javier Escaned

Istituto Cardiovascular, Hospital Universitario, Clinico San Carlos.



- Accurate vessel sizing and length measurement without a pullback device
- Localization of IVUS images with angio co- registration is designed to help reduce the risk of geographic miss, which study data suggests may occur in 66.5% of PCIs¹
- Use with fast, plug and play digital IVUS unique to Philips IGT-D



- Precise lesion severity, location and length assessment
- Visualization of pressure gradients to facilitate stent planning
- Use with FFR or the simplified workflow of hyperemia-free iFR

1. Costa et al. Impact of Stent Deployment Procedural Factors on Long-term Effectiveness and Safety of Sirolimus- Eluting Stents (Final results of the Multicenter Prospective STLLR Trial), Am J Cardiol 2008 Jun 15;101(12):1704-11



Clinical data



Guidelines

ESC/EACTS 2018 guidelines ¹ (actual wording)	Class	Level of evidence
IVUS should be considered to assess the severity of unprotected left main lesions	Ila	B
IVUS should be considered to optimize treatment of unprotected left main lesions	Ila	B
IVUS or OCT should be considered in selected patients to optimize stent implantation	Ila	B
IVUS and/or OCT should be considered to detect stent-related mechanical problems leading to restenosis	Ila	C
IVUS or OCT to assess mechanisms of stent failure	Ila	C

Class Ila: “should be considered”,
Class Iib: “may be considered.”

ACC/AHA/SCAI 2011 guidelines ² use IVUS:	Class	Level of evidence
For the assessment of angiographically indeterminate left main CAD	Ila	B
4 to 6 weeks and 1 year after cardiac transplantation to exclude donor CAD, detect rapidly progressive cardiac allograft vasculopathy, and provide prognostic information	Ila	B
To determine the mechanism of stent restenosis	Ila	C
For the assessment of non-left main coronary arteries with angiographically intermediate coronary stenoses (50% to 70% diameter stenosis)	Iib	B
For guidance of coronary stent implantation, particularly in cases of left main coronary artery stenting	Iib	B

Class Ila: “is reasonable”,
Class Iib: “may be considered.”

1. Neumann et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. EuroIntervention 2019;14:1435-1534


2. Levine G et al, 2011 ACCF/AHA/SCAI Guideline for Percutaneous Coronary Intervention. Circulation 2011;124:e574–e651..

Consensus statements

Society statements advocating the use of IVUS

“For determining optimal stent deployment (complete stent expansion and apposition and lack of edge dissection or other complications after implantation).”

“For determining the size of the vessel undergoing stent implantation”




- Mintz et al. Intravascular ultrasound in the evaluation and treatment of left main coronary artery disease: a consensus statement from the European Bifurcation Club. *EuroIntervention* 2018;14:e467-e474
- Lotfi A, et al. Expert consensus statement on the use of fractional flow reserve, intravascular ultrasound, and optical coherence tomography: a consensus statement of the society of cardiovascular angiography and interventions. *Catheter Cardiovasc Interv*. 2014 Mar 1;83(4):509-18.
- Percutaneous Recanalization of Chronic Total Occlusions: 2019 Consensus Document from the EuroCTO Club. *EuroIntervention* 2019; Jaa-510 2019, doi:10.4244/EIJ-D-18-00826.
- Galassi et al. Utility of Intravascular Ultrasound in Percutaneous Revascularization of Chronic Total Occlusions. *Am Coll Cardiol Intv* 2016;9:1979-91)
- Räber, Johnson et al. Clinical use of intracoronary imaging. Part 1: guidance and optimization of coronary interventions. Part 2: acute coronary syndromes, ambiguous coronary angiography findings, and guiding interventional decision-making: An expert consensus document of the European Association of Percutaneous Cardiovascular Interventions *European Heart Journal* (2018) 0, 1–20. doi:10.1093/eurheartj/ehy285. Part 2. *European Heart Journal* (2019), ehz332, <https://doi.org/10.1093/eurheartj/ehz332>



A preponderance of evidence shows

IVUS benefits patients





Meta-analyses of RCT




Elgendy meta-analysis

IVUS guided DES-implantation was associated with a significant reduction of MACE, cardiovascular mortality and stent thrombosis

- Includes 7 RCTs
- 3,192 patients
- Including only DES patients

At a mean of 15 months, routine IVUS-guided PCI was associated with a significant reduction in:



37% of MACE

(6.5% versus 10.3%; OR, 0.60; 95% CI, 0.46–0.77; P<0.0001)

38% of ischemia-driven target lesion revascularization

(4.1% versus 6.6%; OR, 0.60; 95% CI, 0.43–0.84; P=0.003)

58% of cardiovascular mortality

(0.5% versus 1.2%; OR, 0.46; 95% CI, 0.21–1.00; P=0.05)

54% of stent thrombosis

(0.6% versus 1.3%; OR, 0.49; 95% CI, 0.24–0.99; P=0.04)



Shin meta-analysis

IVUS guided DES-implantation was associated with a significant reduction of MACE

- Meta-analysis with individual patient-level data
- 3 RCTs including high risk coronary anatomy (long lesions or CTO)
- 2,345 randomized patients, all were treated for long lesions or chronic total occlusions (CTO)
- Including only new generation of DES
- Notably, in the present study, the primary endpoint did not include TLR
- IVUS-guided new-generation DES implantation in long lesions and CTO was associated with favourable outcomes for MACE
- MACE was defined as a composite of cardiac death, myocardial infarction, or stent thrombosis





Shin meta-analysis

IVUS guided DES-implantation was associated with a significant reduction of MACE

At 1 year


Major adverse cardiac events had occurred in 0.4% of the patients who underwent IVUS-guided DES implantation versus 1.2% of those who underwent angiography-guided DES implantation.

(Hazard ratio [HR]: 0.36; 95% confidence interval [CI]: 0.13 to 0.99; $p = 0.040$)

MACE by intention-to treat analysis

Hazard ratio: 0.36 (95% CI: 0.13-0.99)

Log Rank: $P = 0.040$





Large randomized control trials



IVUS-XPL

IVUS-guided vs angiography-guided everolimus-eluting stent implantation

- Randomized control trial (RCT)
- 1400 patients treated with EES stent ϵ 28 mm


At 1 year

IVUS-guided group had a 52% reduced risk of MACE compared with the angiography-guided group.

(2.9% vs 5.8%, (hazard ratio [HR], 0.48 [95% CI, 0.28 to 0.83]; P = .007).

These differences were primarily due to lower risk of target lesion revascularization.

All population MACE%





IVUS-XPL


IVUS-guided vs angiography-guided everolimus-eluting stent implantation

Imaging standardization with IVUS defined criteria for optimum stent deployment:

MSA > minimum lumen at distal

Patients who did not meet IVUS criteria for optimum stent deployment had a significantly higher incidence of MACE compared with those who met IVUS criteria (**4.6% vs 1.5%, p = 0.020**)

MACE rate in IVUS-guided group only





IVUS-XPL: 5 year follow up

Compared with angiography-guided stent implantation, IVUS-guided stent implantation resulted in a significantly lower rate of major adverse cardiac events up to 5 years


	IVUS-guidance (n=700)	Angiography-guidance (n=700)	Hazard ratio (95% CI)	Log-rank Pvalue
Primary end point				
MACE	36 (5.6%)	70 (10.7%)	0.50 (0.34–0.75)	0.001
Secondary end point				
Cardiac death	6 (0.9%)	14 (2.2%)	0.43 (0.17–1.12)	0.074
Target lesion related MI	4 (0.6%)	6 (0.9%)	0.67 (0.19–2.36)	0.525
Ischemia-driven TLR	31 (4.8%)	55 (8.4%)	0.54 (0.33–0.89)	0.007
Stent Thrombosis	2 (0.3%)	2 (0.3%)	1.00 (0.14–7.10)	1.000
Acute	1 (0.1%)	1 (0.1%)		
Sub-acute	1 (0.1%)	0		
Late	0	1 (0.1%)		




IVUS-XPL: 5 year follow up

Sustained 5-year clinical benefits resulted from both within 1 year and from 1 to 5 years post-implantation even in the current DES implantation era

At 5 years



Between 1 and 5 years






ULTIMATE trial

Intravascular ultrasound-guided versus angiography-guided implantation of drug-eluting stent in all-comers

- 8 sites, randomized prospective trial, IVUS or angiography guided PCI in China
- 1,448 all-comer patients
- In both arms, all lesions were post dilated with non-compliant balloon inflated at >18 atm

Imaging standardization with 3 IVUS defined criteria, pre & post PCI





- 1** Proximal: Plaque burden at the 5-mm proximal or distal to the stent edge <50%.
- 2** Stent: Expansion satisfactory - minimal lumen CSA in stented segment >5mm², or 90% of distal reference lumen CSA.
- 3** Distal: No edge dissection involving media with length >3mm.

ULTIMATE trial

Intravascular ultrasound-guided versus angiography-guided implantation of drug-eluting stent in all-comers

Primary endpoint based on patient-level comparison



47%

reduced risk of TVF* at one year using IVUS-guided PCI.


*(cardiac death, target vessel MI and clinically driven TVR)



ULTIMATE trial

Intravascular ultrasound-guided versus angiography-guided implantation of drug-eluting stent in all-comers

Primary endpoint for patients who met/didn't meet IVUS criteria



Where optimal IVUS criteria was achieved, patients had the lowest event rates

Only 1.6% TVF at 12 months when optimal IVUS - guided PCI criteria was met in all comer patients



ULTIMATE trial: 3 year follow up

New slide


IVUS guidance associated with a lower risk for 3-year TVF and definite/probable ST

	1Y Follow up				2Y Follow up				3Y Follow up			
	IVUS Guidance (n % 724)	Angiographic Guidance (n % 724)	Hazard Ratio (95% CI)	p Value	IVUS Guidance (n % 724)	Angiographic Guidance (n % 724)	Hazard Ratio (95% CI)	p Value	IVUS Guidance (n % 724)	Angiographic Guidance (n % 724)	Hazard Ratio (95% CI)	p Value
Number of patients	722	722			718	719			714	709		
Target vessel failure	21 (2.9)	39 (5.4)	0.53 (0.31–0.90)	0,02	43 (6.0)	65 (9.0)	0.65 (0.44–0.95)	0,03	47 (6.6)	76 (10.7)	0.60 (0.42–0.87)	0,01
Cardiac death	5 (0.7)	10 (1.4)	0.50 (0.17–1.45)	0,19	9 (1.3)	16 (2.2)	0.56 (0.25–1.26)	0,16	13 (1.8)	19 (2.7)	0.68 (0.34–1.38)	0,28
Target vessel MI	7 (1.0)	11 (1.5)	0.63 (0.25–1.64)	0,34	7 (1.0)	14 (1.9)	0.50 (0.20–1.23)	0,12	7 (1.0)	15 (2.1)	0.46 (0.19–1.14)	0,09
Clinically driven TVR	11 (1.5)	21 (2.9)	0.51 (0.25–1.07)	0,07	31 (4.3)	42 (5.8)	0.72 (0.45–1.15)	0,17	32 (4.5)	49 (6.9)	0.64 (0.41–1.00)	0,05
All-cause death	10 (1.4)	17 (2.3)	0.58 (0.27–1.28)	0,17	24 (3.3)	27 (3.8)	0.88 (0.51–1.53)	0,65	31 (4.3)	31 (4.4)	0.99 (0.60–1.63)	0,98
Clinically driven TLR	9 (1.2)	19 (2.6)	0.47 (0.21–1.03)	0,05	26 (3.6)	40 (5.6)	0.63 (0.39–1.04)	0,07	27 (3.8)	45 (6.3)	0.59 (0.36–0.94)	0,03
Target lesion failure	20 (2.8)	37 (5.1)	0.53 (0.31–0.92)	0,02	38 (5.3)	63 (8.8)	0.59 (0.39–0.88)	0,01	42 (5.9)	72 (10.2)	0.57 (0.39–0.83)	0,003
Definite/probable ST	1 (0.1)	5 (0.7)	0.20 (0.02–1.70)	0,1	1 (0.1)	7 (1.0)	0.14 (0.02–1.15)	0,03	1 (0.1)	8 (1.1)	0.12 (0.02–0.99)	0,02




ULTIMATE trial: 3 year follow up

IVUS guidance associated with a lower risk for 3-year TVF



IVUS guidance associated with a lower risk for 3-year TVF, particularly for patients with IVUS-defined optimal procedures, relative to angiographic guidance among all comers undergoing second generation DES implantation.







ULTIMATE CKD:

CKD patients** undergoing DES implantation associated with a higher risk of TVF at 12 months: the risk of TVF in the CKD patients can be significantly decreased through IVUS guidance.

**CKD -defined as an estimated glomerular filtration rate (eGFR) $<60 \text{ mL min}^{-1} 1.73 \text{ m}^{-2}$ - was available in 1,443 patients, of whom 723 were in the IVUS guidance group, and 720 were in the angiography guidance group: CKD was present in 349 (24.2%) patients.





Large meta-analysis

Ahn meta-analysis

IVUS-guided PCI was associated with a significant reduction in hard endpoints; MACE & MI

- Included 17 studies covering 26,503 patients
- 12,499 patients had PCI with IVUS guidance
- Comprehensive analysis reflecting DES studies over a decade

IVUS guided DES-implantation was associated with a significant reduction of hard endpoints:

25% reduction of MACE
43% reduction of MI



2005-2014 IVUS in Clinical Research

17 STUDIES (14 OBSERVATIONAL AND 3 RANDOMIZED) INCLUDED IN THE AHN META-ANALYSIS OF OUTCOMES AFTER INTRAVASCULAR ULTRASOUND-GUIDED VERSUS ANGIOGRAPHY-GUIDED DRUG-ELUTING STENT IMPLANTATION



Total IVUS Patients: n=12,499 | Total DES Patients: 26,503 | Studies Included: 17




Large observational studies


ADAPT-DES

Assessment of dual anti-platelet therapy with drug-eluting stents

- Largest study ever conducted with IVUS guidance
- Multi-center global registry with 8583 consecutive patients
- 3349 patients had PCI with IVUS guidance
- 64% Xience / Promus stents



Investigators were asked if and how IVUS changed their procedure?



**Others* category may include a combination of "Higher Pressure", "Under Expansion", "Malapposition", and "Additional Stent".


1. Witzenschnabel B et al. Relationship Between Intravascular Ultrasound Guidance and Clinical Outcomes After Drug-Eluting Stents: The ADAPT-DES Study. *Circulation* 2014;Jan: 129,4:463-470.

2. Maehara et al. Relationship Between Intravascular Ultrasound Guidance and Clinical Outcomes After Drug-Eluting Stents: Two-Year Follow-Up of the ADAPT-DES Study. *Circ Cardiovasc Interv.* 2018;11:e006243. DOI: 10.1161/CIRCINTERVENTIONS.117.006243.

ADAPT-DES

IVUS guided PCI was associated with a significant reduction of cardiac death, and ST in addition to MI and MACE

Relationship between IVUS use and MACE (definite/probable ST, cardiac death, MI) within 2 years



Number at risk:	0	6	12	18	24
IVUS used	3361	3206	3117	2988	1739
No IVUS used	5221	4912	4740	4537	2177

35%

reduction in
MACE at 2 years
(4.9% vs. 7.4%, p<0.001)

Data suggests the benefit of IVUS guidance increases from 1 to 2 years.




ADAPT-DES


IVUS guidance was beneficial in complex and standard* PCI.

*Single vessel, non-left main bifurcation & stable CAD

Association of IVUS use with MACE (definite/probable ST, cardiac death, MI) in relation to lesion complexity




Association of IVUS use with MACE (definite/probable ST, cardiac death, MI) in relation to index presentation




ADAPT-DES

**IVUS guidance was of greatest benefit in complex PCI;
multi-vessel, bifurcation, left main and unstable angina patients**

Association of IVUS use with MACE (definite/probable ST, cardiac death, MI) in relation to lesion complexity




Association of IVUS use with MACE (definite/probable ST, cardiac death, MI) in relation to index presentation





ADAPT-DES IVUS sub study

IVUS guidance reduces cardiac death, and ST in addition to MI and MACE within 2 years after DES implantation



*These calculations have been completed by Volcano and source data from Dr. Witzentichler's 2013 TCT presentation¹.

**MACE defined as definite/probable ST, cardiac death, and MI.

Witzentichler, Bernhard. "Does IVUS Reduce Stent Thrombosis with DES? Two-year results from the prospective, multicenter ADAPT-DES study." TCT 2013. Moscone Center, San Francisco, CA. 29 October 2013.



DEFINITION registry

IVUS-guided drug-eluting stent implantation is associated with improved clinical outcomes in patients with unstable angina and complex coronary artery true bifurcation lesions

- Prospective, multi-center registry
- 1465 patients with unstable angina and Medina 1,1,1 & 0,1,1 coronary bifurcation
- Clinical follow up period of up to 7 years
- Importance to reach standardized IVUS criteria to obtain optimal clinical benefit.

Standardized IVUS criteria:

- Stent symmetry index > 70%
- Stent expansion > 90%
- No malapposition,
- No Type B-Cdissection
- No impairment of bloodflow

Patients without optimal IVUS-guided stent expansion had worse outcome than those with optimal IVUS-guided results.


MACE 16.9% vs 8.6% (p=0.029)



DEFINITION registry

Reported evidence showing that IVUS benefits patients with complex coronary artery true bifurcation lesions


At 1 year



IVUS-guidance was associated with a significant reduction in Major Adverse Cardiac Events (MACE)

MACE in IVUS group occurred in 10.0% vs 15.0% (p = 0.036)

At 7 years



IVUS-guidance resulted in lower:



- **Cardiac death rate (6.5 versus 1.3%, p = 0.002)**
- **MI (8.4 versus 2.3%, P < 0.001)**



SCAAR left main registry

IVUS-guidance is associated with better outcome in patients undergoing unprotected left main coronary artery stenting compared with angiography guidance alone

- Prospective, Swedish nationwide observational study.
- 2468 patients who underwent unprotected LMCA PCI between 2005 and 2014 due to CAD or ACS.
- Clinical follow up period of up to 10 years
- Importance to reach appropriate stent size as larger stents were independently associated with improved outcome.




The improvement in outcome was primarily driven by a marked and significant reduction of all-cause mortality (-35% death in total population, -46% death in 340 propensity pair matched population).



Comparison of left main observational studies

IVUS guidance associated with reduced mortality and MACE rates in LMCA stent implantation^{1,2,3}

All courses of mortality IVUS in left main



Reduction of all cause of mortality consistently found in all observational studies

1. Park SJ, et al. MAIN COMPARE. Circ Cardiovasc Interv. 2009;2(3):167-77. 2. de la Torre Hernandez JM, et al. IVUS-TRONCO-ICP, JACC Cardiovasc. Interv. 2014;7(3):244-54. Andell et al. IVUS in LMCA PCI, Circ Cardiovasc Interv 2017; 10(5).

A preponderance of evidence shows that IVUS benefits patients

IVUS benefits patients

Using a standardized IVUS criteria supports optimal clinical benefit

IVUS-XPL, Hong et al. 2015 ULTIMATE, Zhang J et al. 2018 DEFINITION registry, Chen et al. 2017

Multiple studies show patients without optimal IVUS-guided stent expansion had a significantly worse outcome than those with optimal IVUS-guided results.


DEFINITION registry, Chen et al. 2017
MACE 16.9% vs 8.6% (p=0.029)

Defined criteria for ULTIMATE Zhang J et al. 2018

1. Plaque burden at the 5-mm proximal or distal to the stent edge <50%.
2. Expansion Satisfactory - Minimal lumen CSA in stented segment >5mm², or 90% of distal reference lumen CSA.
3. No edge dissection involving media with length >3mm.

Only 1.6% TVF at 12 months when optimal IVUS-guided PCI criteria was met in all-comer patients.

Primary endpoint for patients who met/didn't meet IVUS criteria





Glossary

ACC	American College of Cardiology
AHA	American Heart Association
ATM	Atmosphere
CAD	Coronary Artery Disease
CSA	Cross-Sectional Area
DES	Drug Eluting Stent
EACTS	European Association for Cardio-Thoracic Surgery
EAPCI	European Association of Percutaneous Cardiovascular Interventions
EBC	European Bifurcation Club
EEM	External Elastic Membrane
EP	Electrophysiology
ESC	European Society of Cardiology
FFR	Fractional Flow Reserve
GM	Geographical Miss
ICE	Intracardiac Echocardiography
iFR	Instantaneous wave-Free Ratio
IGT	Image Guided Therapy
IGT-D	Image Guided Therapy – Devices
IGT-S	Image Guided Therapy – Systems

ISR	In-Stent Restenosis
IVUS	Intravascular Ultrasound
LMCA	Left Main Coronary Artery
MACE	Major Adverse Cardiovascular Event
MHz	MegaHertz
MI	Myocardial Infarction
MLA	Minimum Lumen Area
MLD	Minimum Lumen Diameter
OCT	Optical Coherence Tomography
PCI	Percutaneous Coronary Intervention
PV	Peripheral Vascular
RCA	Right Coronary Artery
RCT	Randomized Control Trial
RLA	Reference Lumen Area
SCAI	Society for Cardiovascular Angiography and Interventions
ST	Stent Thrombosis
TLR	Target Lesion Revascularization
TVF	Target Vessel Failure
TVR	Target Vessel Revascularization

