

PHILIPS

MX750

Mid-acute care



**Supporting productivity
and flexibility**
in mid-acute care settings

A monitoring solution that meets your clinical, IT and patient needs

Patients in mid-acute care have different monitoring requirements to critical and anesthesia cases, and the staff looking after them have different needs and priorities too. The Philips IntelliVue MX750 for mid-acute care has been specially adapted for these care settings. It offers exceptional functionality, as well as a streamlined, flexible workflow and dedicated, scalable IT and clinical services that help you unlock the full value and potential of your equipment.

Simple to use, easy to see

The MX750 for mid-acute care has a full HD, 19" projected capacitive (pCAP) touchscreen for ease of use and excellent readability. The large display means that physiological data is visible from a distance, helping staff identify changes in condition even when they are not directly at the bedside.

The device offers a bed-to-bed overview, which allows caregivers to respond to alarms swiftly – even if they are currently tending to another patient. Particularly in settings where there are fewer staff per patient, this helps support early intervention in the event of deterioration. Optionally, the IntelliVue MX750 for mid-acute care offers Early Warning Scorecards and a highly configurable EWS score generator.

Compact and secure

The monitor is a compact, single-device solution with one power plug. It is highly resistant to disinfection agents. Thanks to simple barcode scanning via a standard USB scanner, the IntelliVue MX750 for mid-acute care offers easy patient admission and integration into hospital information systems.

The solution has been developed based on the 'security by design' philosophy, which proactively addresses risks and threats, making it a robust addition to your IT infrastructure. As a mid-acute solution, the IntelliVue MX750 offers a scalable software-based upgrade path, should you need to add high-end capabilities or deploy your devices in more critical-care environments in the future.



Simplifying clinical workflow

Philips IntelliVue MX750 patient monitors for mid-acute care conveniently bridge the clinical and IT worlds to simplify workflow for clinicians, by giving you integrated access to the solutions that help you deliver efficient and effective care.

With smooth surfaces, minimal seams, and tested for use with a variety of approved cleaning agents, IntelliVue MX750 monitors for mid-acute care are simple to clean. The device's fan-less operation and passive cooling support your efforts to reduce cross-contamination.

A familiar interface


The IntelliVue MX750 for mid-acute care features an intuitive user interface that is consistent with your existing IntelliVue patient monitors. This helps reduce training effort and supports staff productivity.

Make the most of your EMR


The IntelliVue MX750 for mid-acute care lets you enhance your investment by helping reduce cost and complexity when connecting bedside devices to your EMR. If you are using the optional IntelliVue XDS software, you can even view the patient monitor remotely, while accessing the patient's record simultaneously on the EMR. Moreover, these patient monitors offer the optional capability to display native Web applications and Citrix®-published applications.

Extensive integration support


Alongside these technological features, we can also support you in integrating our solutions into your network. We can draw upon over 30 years of experience and expertise in the field of healthcare IT – meaning that whatever your integration challenges, we are ready and able to help you meet your needs. Moreover, we offer a comprehensive range of clinical services to help you extract the maximum from your technology by building up the knowledge, education and skills of your care delivery teams.



MX750 for mid-acute care provides a convenient overview of multiple beds.



The device gives a clear overview of numeric trends.



The device provides guidance on correct ECG lead placement.


Robust cybersecurity

The IntelliVue MX750 patient monitor for mid-acute care offers a host of cybersecurity capabilities, including certificate-based node authentication, plus the FIPS 140-2 certified encryption of network data, file systems, print reports, remote display data and remote control commands. In short, the monitor gives you the same level of protection that the IntelliVue MX750 for anesthesia and critical care delivers.



Be alert to subtle yet meaningful clinical changes


The IntelliVue MX750 for mid-acute care offers a host of clinical features that let you concentrate on relevant changes in your patient’s condition. Optionally, the device can deliver the full range of clinical features included in the high-end version of the IntelliVue MX750.



The 19" display provides clarity and an intuitive graphical representation of changes in vital signs.




(Optional) The Alarm Advisor application encourages critical thinking for deploying alarm management strategies or personalizing alarm thresholds depending on the patient’s condition (according to hospital-defined guidance).




(Optional) Early Warning Score (EWS) helps you identify relevant changes in patient condition early on.



Creating event-trigger groups, using Advanced Event Surveillance, to get higher priority alerts when multi-parameter events occur.



(Optional) Continuous ST segment monitoring may identify early warning signs of coronary ischemia.¹



The pop-up bed-to-bed overview alerts you to alarms on another patient’s monitor.



Support and services

Pairing equipment with the right services helps you get more from your systems across their lifecycle. Philips Services encompass everything from service contracts that help you keep your installation at the cutting edge, to lifecycle management and maintenance services IT and project management services that allow you to extract the full value of your equipment throughout the lifecycle. We also offer clinical services that empower you to optimize your care environment, and IT services, that may support you with your interoperability challenges.

Protecting investments, protecting your systems

The IntelliVue MX750 for mid-acute care is a next-generation platform that allows you to utilize a wide range of features and functions through software and hardware maintenance agreements that can be tailored to fit your specific needs. Philips can help you maintain your solutions and keep them up-to-date and secure. For example, we can help you bring the latest software versions to your systems – supporting you in equipping your organization in line with the latest security standards, while helping you safeguard against potential vulnerabilities.

At the same time, we also understand how important it is to keep your staff up-to-date on the latest clinical and IT developments, and fully aware of new functionalities. That's why we offer a range of professional services that support you in unlocking the true value of your investments, by helping your staff to build the right knowledge, education and skills.



To find out more about the Philips IntelliVue MX750 for mid-acute care, please contact your Philips representative.

Reference

1. Sandau et al. Update to Practice Standards for Electrocardiographic Monitoring in Hospital Settings: A Scientific Statement From the American Heart Association. *Circulation*.2017;136:e273-e344.doi:10.1161/CIR.0000000000000527

