

PHILIPS

sonicare



Cordless water flosser
with accessories

Power Flosser 3000

Quad Stream technology

Pulse waves guide as you go



HX3806/31

Flossing reinvented with QuadStream technology

Removes up to 99.9% of plaque* in treated areas

Effortless thorough cleaning between teeth to improve your gum health. Unique Quad Stream tip cleans more surface area with less effort for more effective flossing every time, so it's easy to get flossing right.

Easy and convenient

- Easy-to-fill 250 ml reservoir
- Effortless technique
- Standard tip for spot cleaning
- Fast, universal charging with USB-C cable

Faster, more effective flossing

- 2 flossing modes, 3 intensities
- Complete clean in 60 seconds
- Pulse Wave technology guides from tooth to tooth
- Quad Stream technology for faster, more effective flossing

Power Flosser 3000

Cordless water flosser with accessories

HX3806/31


Highlights

Quad Stream technology



Unique X-shaped Quad Stream tip separates flow into 4 water jets, which cover more surface area between teeth and along the gum line. Get a faster, deeper clean without the hassle of string floss.

Pulse Wave technology



Gentle pulses of water keep you going and guide you from tooth to tooth in Deep Clean mode, so you get it right every time.

2 modes, 3 intensities



Choose the clean that suits your smile. Clean mode delivers continuous flow for effective everyday cleaning. Deep Clean pulses for more thorough cleaning. Adjustable intensity level adapts for comfort.

Complete clean in 60 seconds



Access those hard-to-reach areas with a nozzle that rotates 360 degrees, so that you can use the device in any orientation. Get a thorough 360° clean in just 60 seconds from start to finish.

USB-C Charging



Convenient charging with a USB-C cable uses the same connection as most mobile phones or tablets. Get up to 14 days of cleaning from a single charge.

Standard tip for spot cleaning



A standard single-stream tip boosts pressure for spot cleaning and removing food debris when you need it.

Specifications

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Nozzle attachment: Easily clicks on and off

Performance

Cleaning: Clean entire mouth in 60–90 sec

Items included

Cordless Power Flosser: 1

F1 Standard nozzle: 1

F3 Quad Stream nozzle: 1

Travel bag: 1

USB charging cable: 1

USB wall adapter: 1

Modes

Clean: For exceptional everyday cleaning

Deep Clean+: For an invigorating deep clean

Intensities: 3

Health benefits

Plaque removal: Removes up to 99% of plaque in treated areas*

Design and finishing

Colour: White

