

PHILIPS

sonicare



Interdental -
Rechargeable

AirFloss

w/ 1 nozzle



HX8221/02

Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

Improves oral health

- Helps improve gum health in just two weeks
- Helps prevent cavities between teeth
- Removes up to 5x more plaque between teeth*

Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

Easy to use


- Technology that cleans with a touch
- Guidance tip ensures correct placement
- Comfortable ergonomic handle

Helps build a healthy habit

- An easy way to improve cleaning between teeth

Highlights

Air and micro-droplet technology




AirFloss utilizes unique air and micro-droplet technology to send a burst of pressurized air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

30 Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

Easy-fill reservoir




Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

Guidance tip




The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

Ergonomic handle



Designed to be easy to maneuver in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.

Helps prevent cavities




By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

Removes more plaque



Cleaning between your teeth is an important part of your oral health. Sonicare AirFloss uses air and micro-droplet technology to remove up to 5x more plaque between teeth than manual brushing alone.

An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Specifications

Power

Voltage: Multi-voltage charger

Technical specifications

Battery: Rechargeable

Operating time: Replace nozzle every 6 months

Battery type: Lithium ION

Charge time: 24 hours to full charge

Design and finishing

Color: Pink

Nozzle color: Pink with white tip

Ease of use

Battery indicator: Flashing indicates low battery

Handle: Slim ergonomic design

Nozzle attachment: Easily snaps on and off

Battery Life: Up to 3 wks when fully charged

Items included

Charger: 1 multi-voltage charger

AirFloss Nozzle: 1

Nozzle holder: 1

Cleaning performance

Health benefits: Improves gum health in only two weeks

Performance: Removes up to 5x more plaque*

Nozzles: Best results change every 6mos

Auto-burst feature: Cleans entire mouth in 30 sec

