

PHILIPS



InfraCare infrared
lamp

300W

Zone treatment area

Digital timer



HP3631

Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 300 Watt infrared lamp is very comfortable and the warmth penetrates deep into the skin, stimulates blood circulation and gives relief for "zone" areas of 40x30 cm.

Effective pain relief

- Deep penetrating Infrared warmth

Highlights

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created

by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Specifications

Effective pain relief

Deep penetrating warmth

Zone treatment

Treatment area: 40x30 cm

Infrared halogen lamp: 300 W

Treatment control

Digital timer: with auditative feedback

Technical specifications

Power: 300 W

Battery: LR-54-Alkaline, 1,5 V (lifetime 200 actions)

Cord length: 200 m

Voltage: 230 V

Frequency: 50/60 Hz

Insulation: Class II (double isolation)

Lifetime of lamps: 500 hour(s)

Weight and dimensions

Product weight: 1.4 kg

Product dimensions: 20x30.9x17.2 (WxHxD) cm

A-box dimensions: 29x40.4x74.5 (WxHxD) cm

F-box weight: 1.8 kg

No. of F-boxes in A-box: 3 pcs

Qnt. on Euro pallet: 60 pcs

A-box weight: 5.4 kg

F-box dimensions: 28x38.6x24.5 (WxHxD) cm

Logistic data

Country of origin: Germany

CTV code: 884363101000

Safety

IEC certified: Complies to IEC 60601 and 60335

Easy to use

Adjustability: 0-40 degrees (backwards)

On/off switch

Soft-touch handgrip: for easy transport and positioning

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

