

PHILIPS

sonicare



Sonic electric
toothbrush

2100 Series

Sonic technology

QuadPacer and SmarTimer

Slim ergonomic design

14-day battery life



HX3651/12

Goodbye manual toothbrush. Hello Sonic technology.

Up to 3x better plaque removal*

Sonic technology combined with our brush action gently removes plaque up to 3x better* than a manual toothbrush.

Designed around you

- Ergonomic design makes the toothbrush easy to hold and use
- 14-day long lasting battery life

Optimize your brushing

- Optimize your brushing with SmarTimer and QuadPacer

Easy transition to Sonicare

- Easy-Start to ease the transition

Protects your gums

- Safe and gentle to use

Provides superior clean

- Up to 3x better plaque removal vs a manual toothbrush*
- Our unique technology brings you a powerful yet gentle clean


2100 Series

Sonic electric toothbrush

HX3651/12

Highlights

Advanced Sonicare Technology



Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two month's worth of manual brushing in just 2 minutes.** 31000 brush strokes per minute gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

Removes plaque 3x better*



The Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove plaque up to 3X better* than a manual toothbrush. It removes plaque from your teeth and along the gum line while protecting your gums.

Safe and gentle to use



Sonicare toothbrush provides you with a thorough clean without being harsh on your teeth and gums. Gentle yet powerful sonic vibrations provide an exceptional clean while being gentle on your teeth and gums.

SmartTimer and QuadPacer



2 minute SmarTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

Designed around you




Slim, light-weight, and ergonomic design makes the toothbrush easy to hold and use so that you can get an amazing clean without any strain.

Long lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

Easy-Start



Our Easy-Start program gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

2100 Series

Sonic electric toothbrush

HX3651/12

Specifications

Power

Voltage: DC5V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): 14 days*****

Battery type: Lithium ION

Design and finishing

Color: Light blue

Service

Warranty: 2-year limited warranty

Ease of use

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Items included

Charger: 1

Handle: 1 2100 Series

Brush head: 1 C1 ProResults standard

Cleaning performance

Performance: Removes up to 3x more plaque*

Speed: Up to 31,000 brush strokes per minute

Timer: Quadpacer and SmarTimer

Modes

Clean: For exceptional everyday clean

Consumer Trade Item

Height: 22 cm

Width: 10.9 cm

Length: 6.3 cm

Net weight: 0.184 kg

Gross weight: 0.264 kg

GTIN: 8710103985495

Country of origin: CN

Harmonized system code: 85098000

Outer Carton

Height: 25.4 cm

Width: 20.8 cm

Length: 12.5 cm

Net weight: 0.552 kg

Gross weight: 0.946 kg

GTIN: 18710103985492

