

PHILIPS


AVENT



Comfort Double electric breast pump

Two 4oz bottles and nipples
+travel bag

More milk in less time



SCF334/04

More comfort, more milk

Breast pump with soft massage cushion

When you are comfortable and relaxed, your milk flows more easily. That is why we created our most comfortable breast pump yet: Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow.

Sit more comfortably with no need to lean forward

- More comfortable pumping position due to unique design

Gently stimulates let down and milk flow

- Soft massage cushion with warm feel

Other benefits

- Includes our Natural bottle and nipple for natural latch on
- Intuitive assembly. Easy visual matching of parts
- Compact lightweight design, includes handy travel bag
- Easy cleaning due to the small number of separate parts
- Bonus: Breast pads: 20 day, 20 night, travel pouch, bra clip

Choose your most comfortable setting

- Features a gentle stimulation mode and 3 pumping settings

More milk in less time

- The perfect time saver for moms

Highlights

More comfortable position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle, even when you are sitting up straight. This means you can sit more comfortably when pumping: no need for you to lean forward to make sure all your milk ends up in the bottle. Sitting comfortably and being relaxed when pumping, naturally helps your milk to flow more easily.

Simple settings to choose from

When switched on, the pump automatically starts in gentle stimulation mode to get your milk flowing. Then choose from 3 pumping settings to make milk flow most comfortable for you.

Soft massage cushion



Our massage cushion has a new soft velvety texture that gives a warm feel to the skin for comfortable, gentle stimulation of your milk flow. The cushion is designed to gently mimic your baby's suckling to help stimulate let down.

Easy cleaning

Cleaning is easy, thanks to the small number of separate parts. Your milk will never come into contact with the tubing and base unit. All parts are dishwasher proof, except the electrical parts.

Perfect time saver for moms

Simultaneous pumping at both breasts is proven more efficient and may even boost your ability to produce breast milk.*

Intuitive assembly

Easy visual matching of parts for intuitive assembly.

Compact design

The breast pump has a compact design, which makes it easy to hold and position on your breast. The small, lightweight base unit can easily be placed within comfortable reach for full control when pumping. For extra convenience in transport and storage, the tube simply wraps around the base unit. The breast pump comes with a handy travel bag.

Bonus items

Includes 20 day breast pads, 20 night breast pads, 1 travel pouch and 1 breastfeeding reminder (bra clip).

Includes our Natural nipple

The wide breast shaped nipple promotes natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Specifications

Breast pump

Material: BPA free*

Design

Bottle design: Ergonomic shape, Wide neck
Breast pump design: Compact design

Country of origin

England

Material

Breast pump: Polypropylene, BPA free*

Bottle: Polypropylene, BPA free*

Nipple: Silicone, BPA free*

What is included

Breast pad sample packs: 2 (2 Day pads & 2 Night pads) pcs

Extra soft newborn flow nipple: 2 pcs

Sealing disc for milk storage: 2 pcs

Natural bottle 4oz: 2 pcs

Breast pump body: 2 pcs

Standard size cushion: 2 pcs

Base unit incl. tubing: 1 pcs

Travel bag: 1 pcs

Spare diaphragm: 2 pcs

Travel cover: 2 pcs

Bonus day breast pads: 20 pcs

Bonus night breast pads: 20 pcs

Bonus travel pouch: 1 pcs

Bonus bra clip: 1 pcs

Ease of use

Breast pump use: Intuitive assembly, Easy cleaning, Fully compatible range

Functions

No leaning forward: Sit in a comfortable position

Soft massage cushion: Gentle stimulation

Settings: 3 Expression settings, 1 Stimulation mode

More milk in less time: The perfect time saver

Development stages

Stage: 0 - 6 months

